





**Become part of the biofilm  
success story**




*A new technology*

[contact@biofilmcontrol.com](mailto:contact@biofilmcontrol.com)

# The **BioFilm Ring Test**<sup>®</sup> is a fast and convenient new approach for studying biofilm formation

This unique and proprietary approach monitors how paramagnetic beads are immobilized during the formation of the biofilm:

- ⇒ A magnet is used to assemble the remaining non-immobilized beads into a macroscopic spot.
- ⇒ The resulting spot is quantified through specialized image algorithms.




**This technology is protected worldwide by several patents.**

The **BioFilm Ring Test**<sup>®</sup> is operated through the **BioFilm Scan M** using ready-to-use kits provided exclusively by **BioFilm Control**.

The **BioFilm Scan M** includes:

- A precision 96 well magnetic plate holder,
- A specialized scanner,
- A ready-to-use high quality desktop computer,
- The data-acquisition and analysis software already installed and ready to run.



Each kit contains

- 10 microtiter plates or 120 strip-wells in SBS format specially designed and manufactured for biofilm applications.
- Negatively or positively charged paramagnetic beads.
- A special anti-reflective solution used at image acquisition time to allow reliable image quantification.

## The procedure is extremely simple and fast to operate

- a) Add the paramagnetic beads to the bacterial culture solution, mix and dispense it into the wells
- b) Incubate the desired time at the desired temperature.
- c) Add the anti-reflective solution
- d) Make a first reading using the **BioFilm Scan M**,
- e) Place the wells on the magnetic holder
- f) Make a second reading using the **BioFilm Scan M**,
- g) Let the software analyze the images and quantify the extend of immobilized beads
- h) Proceed further with data analysis: data is provided directly in a standard spreadsheet file

The **BioFilm Ring Test**<sup>®</sup> specific time overhead is of the order of 5 minutes, even for a full 96 well micro-titer plate.




Fig. 1: the **BioFilm Ring Test**<sup>®</sup> workflow

## The **BioFilm Ring Test**<sup>®</sup> has been thoroughly validated

The proof of concept has been conducted in collaboration with a group of INRA (the French national research institute for agronomic research) and have been published in a peer-reviewed journal: Chavant *et al.*, Journal of Microbiological Methods (2007) **68**, 605-12.

For instance, scanning electron microscopy clearly shows that the beads are immobilized by the biofilm matrix of *Listeria monocytogenes* (Fig. 2.A) and *Staphylococcus xylosus* (Fig. 2B), which are respectively low and high EPS producers.





Fig. 2: scanning electron microscopy images(false colors)showing paramagnetic beads (red-brown spheres) immobilized on *Listeria monocytogenes* (grey sticks, A) and *Staphylococcus xylosus* (yellow spheres, B) biofilms.

The technology has already been applied on a wide variety of **bacterial and algal micro-organisms**:

- ✓ *Escherichia coli*
- ✓ *Staphylococcus*
- ✓ *Pseudomonas*
- ✓ *Lactobacillus*
- ✓ *Salmonella*
- ✓ *Streptococcus*
- ✓ *Arthrospira*
- ✓ *Rhodella*
- ✓ *Listeria*
- ✓ ...

## The key strengths of the **BioFilm Ring Test**<sup>®</sup>

The Cristal Violet has been for many years the reference method used for studying biofilm formation. Compared to this reference method (Fig. : 3), the **BioFilm Ring Test**<sup>®</sup>:

- ⇒ Gives the **same qualitative results**
- ⇒ Is **significantly more sensitive**: bacterial adhesion and the initial steps of **biofilm formation can be detected and scored already after 15 minutes incubation**
- ⇒ Requires **much less manipulations**
- ⇒ Can be **fully automated using high-throughput robots**
- ⇒ Is **much faster to operate, overhead is of the order of 5 minutes**, even for a full 96 well plate
- ⇒ **Use only reactive safe for the user**
- ⇒ Gives information about the **early stage of biofilm formation, including surface conditioning and primo adhesion.**




Fig. 3 : biofilm formation of *E.c* DH5α (violet), *S.c* TM300 (red) *S.x* DSM 20267 (blue) and *L.m* EGDe (green) monitored using the Cristal Violet method (top) and the **BioFilm Ring Test**<sup>®</sup> (bottom).

## The **BioFilm Ring Test**<sup>®</sup> has already been used for a wide variety of applications:

- ⇒ High-throughput screening of **curative and preventive antibacterial agents and antibiotics**
- ⇒ Characterization of **anti-biofilm coating agents**
- ⇒ Research of **pathogenicity/virulence signatures** of micro-organisms
- ⇒ Study of the **effect of enzymatic activity on biofilm formation and degradation**, for instance the action of dextranases, **proteases** and **DNAses**
- ⇒ Study of the **degradation of purified exo-polysaccharides by enzymes**
- ⇒ Study of **molecules secreted by micro-organisms** (when such molecules induce the immobilization of the paramagnetic beads to the surface of the well)
- ⇒ Measuring the **adherence of the biofilm** on the surface...

## The **BioFilm Ring Test**<sup>®</sup> is also available as a service


- ⇒ We offer **several years of expertise and innovation** in biofilm analysis,
- ⇒ We operate routinely in a **class P2 laboratory**, and can access to a partner's **P3 class laboratory** when required
- ⇒ We can perform **complementary analyses** such as Cristal Violet, fluorescence microscopy, ...
- ⇒ We have developed a **unique high-throughput robotic system** performing the **BioFilm Ring Test**<sup>®</sup> in a fully automated manner on large number of samples and numerous conditions.
- ⇒ BioFilm control is on the track of becoming **ISO9001:2008 certified** in January 2010.




## The **BioFilm Ring Test**<sup>®</sup> is leading to success

This method is used by a growing number of industrial and academic groups. It has generated already 7 peer-reviewed articles (see attached list for full details) and numerous posters presented at international conferences (full details and copies on request):


 *Journal of Applied Microbiology* (2009)

 *J Mol Microbiol Biotechnol.* (2009)

 *PLoS One.* (2008)

 *Int J Antimicrob Agents.* (2009)

 *Invited presentation at CHRO2009* (2009, Japan)

 *Posters at: Eurobiofilm2009* (2009, Italy), *2nd Meeting of IADR Pan Asian Pacific Federation* (2009, China), *4th IDF Dairy Science & Technology Week*, (2009, France), *3rd International Conference Biofilm III* (2008, Germany), ...

 *Int J Food Microbiol.* (2009)

 *Appl Biochem Biotechnol.* (2008)

 *J Microbiol Methods.* (2007)



**Become part of the biofilm success story...Contact us:**

**contact@biofilmcontrol.com**

**+33 (0) 473 333 980**

---

### BioFilm Control

SAS au capital de 790.500 € - SIRET 480 060 284 00019 – RCS RIOM – TVA intracommunautaire FR39480060284  
Biopôle Clermont-Limagne, 63360 Saint-Beauzire, FRANCE  
Tel : +33 (0) 473 333 380 - Mél : contact@biofilmcontrol.com