

POST-OPERATIVE SHOE WITHOUT FOREFOOT SUPPORT "BAROUK" and "MAYZAUD" long sole version

DESCRIPTION

post-opérative shoe without forefoot support

is available in two versions : ordinary version and a new version with a foreplate

The **ordinary version** consists in :

- an upper, made of a strong navy-blue jean strengthened in the heel part by a special stiffener for guaranteeing the seating of the backfoot.
- three velcros straps serving as fasteners
- heel made of an anti-shock and compressible material including a shock absorber in the front part and offering a 10° talus incline

In addition, the **new version** offers a foreplate which protects the forefoot and avoids the medio-tarsal pain while holding forefoot relief and an inner surface completely padded and without seams

The sizes of the foreplate and its materials guarantee a normal and correct step rolling

INDICATIONS

- principally immediately after any type of forefoot surgery
- all forefoot problems forbidding any forefoot support (inflammation, infection, pain, oedema of the forefoot, injuries or fractures of toes or metatarsal heads or necks)
- algodystrophic syndrome
- diabetic foot ulcers
- rheumatic or dystrophied foot

FUNCTION

This orthosies relieves the forefoot by avoiding any floor contact of the anterior aspect of the foot through a 10° dosiflexed inclination of the heel.

CONTRA-INDICATIONS

- after any type of heel or middle-foot surgery
- all calcaneum problems

SIZING INFORMATION

Each version of the post-operative shoes without forefoot support is available in 4 sizes.:

Size	Pariser Stick Ordinary version	Pariser Stick New version with a foreplate
1	36 - 39	36 - 38
2	39.5 - 42.5	39 - 42
3	43 - 46	43 - 46