

HEEL RELIEF SHOE " SANITAL "

DESCRIPTION

- 10° Plantarflexed inclination of the weight
- bearing sole
- sole : elastomer
- internal sole : polyethylene and latex foam with rounded transition edges
- weight-bearing sole begins just after the heel area
- velcros straps – possibility of foot lifters
- navy blue jean

INDICATIONS

- Any medical pathological condition requiring weight relief on the heel area
- heel ulcers, escharoses
- heel burns, troubles
- heel decubitus
- algodystrophies
- heel spurs
- rheumatic troubles, all the heel problems
- fractures of the talocalcaneal joint

FUNCTIONS

Worn on one foot or on the both feet simultaneously, the Sanital shoe provides new support points and modification of the spreading of pressures.

- complete heel relief by avoiding any floor contact of the posterior aspect of the foot through a 10° plantarflexed inclination and a shortening of the weight bearing sole
- well seating of the ankle in order to reinforce the proprioceptive reflexes
- normal step rolling without restricting the foot function

- early patient mobility : with the use of this shoe long term confinement to bed is prevented, thus hospitalisation time is reduced and the healing process is accelerated

- if worn on one foot, the low height of the outsole makes the walking easy and avoid to limp.

ADVICE FOR USE

Sanital shoe can be worn on one foot or on both feet simultaneously. The foot must be put on in a sock or a protecting jersey in order to avoid skin chafing

SIZE INFORMATION

Each version of Sanital Heel Relief Shoe is available in 4 sizes

Size	Pariser Stitch
1	35/37
2	38/40
3	41/43
4	44/46